



Human JMJD8 peptide (DAG-P1880)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Sequence Similarities	Contains 1 JmjC domain.
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	JMJD8 jumonji domain containing 8 [Homo sapiens (human)]
Official Symbol	JMJD8
Synonyms	JMJD8; jumonji domain containing 8; PP14397; C16orf20; jmjC domain-containing protein 8; jumonji domain-containing protein 8;
Entrez Gene ID	339123
mRNA Refseq	NM_001005920.2
Protein Refseq	NP_001005920.2
UniProt ID	Q96S16
Chromosome Location	16p13.3