



Human TSPAN33 peptide (DAG-P1913)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Plays an important role in normal erythropoiesis. It has a role in the differentiation of erythroid progenitors.
Specificity	Predominantly expressed in erythroblasts.
Purity	70 - 90% by HPLC.
Conjugate	Unconjugated
Sequence Similarities	Belongs to the tetraspanin (TM4SF) family.
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	TSPAN33 tetraspanin 33 [Homo sapiens (human)]
Official Symbol	TSPAN33
Synonyms	TSPAN33; tetraspanin 33; PEN; PEN.; TSPAN-33; tetraspanin-33; penumbra; proerythroblast new membrane;
Entrez Gene ID	340348
mRNA Refseq	NM_178562.4
Protein Refseq	NP_848657.1

UniProt ID	Q86UF1
-------------------	--------

Chromosome Location	7q32.1
----------------------------	--------