



Human SUPT6H peptide (DAG-P1933)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Acts to stimulate transcriptional elongation by RNA polymerase II.
Specificity	Ubiquitously expressed.
Conjugate	Unconjugated
Sequence Similarities	Belongs to the SPT6 family.Contains 1 S1 motif domain.Contains 1 SH2 domain.
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	SUPT6H suppressor of Ty 6 homolog (S. cerevisiae) [Homo sapiens (human)]
Official Symbol	SUPT6H
Synonyms	SUPT6H; suppressor of Ty 6 homolog (S. cerevisiae); SPT6; SPT6H; emb-5; transcription elongation factor SPT6; hSPT6; tat-CT2 protein; tat-cotransactivator 2 protein; histone chaperone suppressor of Ty6;
Entrez Gene ID	6830
mRNA Refseq	NM_003170.3
Protein Refseq	NP_003161.2
UniProt ID	Q7KZ85
Chromosome Location	17q11.2

Function

DNA binding; histone binding; hydrolase activity, acting on ester bonds; poly(A) RNA binding; protein binding; sequence-specific DNA binding transcription factor activity;
