



IAV HA blocking peptide (DAG-P2054)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	BL
Format	Liquid
Buffer	Preservative: 0.02% Sodium Azide Constituents: 0.1% BSA, PBS, pH 7.2
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Preservative: 0.02% Sodium Azide Constituents: 0.1% BSA, PBS, pH 7.2

BACKGROUND

Introduction	Influenza A virus causes influenza in birds and some mammals, and is the only species of influenza virus A. Influenza virus A is a genus of the Orthomyxoviridae family of viruses. Strains of all subtypes of influenza A virus have been isolated from wild b
Keywords	Hemagglutinin; HA; haemagglutinin; IAV HA; Influenza A Hemagglutinin; ha; HA1; HA2; Hemagglutinin ha1 chain; Hemagglutinin ha2 chain; Hemagglutinin precursor