



His tag peptide (DAG-P2055)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

BACKGROUND

Introduction	The H-H-H-H-H-H motif is used as a tag on many recombinant proteins to facilitate purification. The antibody recognizes the His-tag fused to the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.
Keywords	6 His tag; 6 His epitope tag; Hexa His tag; HHHHHH epitope tag; HHHHHH tag; His tag; Polyhistidine Tag