



LAMP1 peptide (DAG-P2064)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Presents carbohydrate ligands to selectins. Also implicated in tumor cell metastasis.
Conjugate	Unconjugated
Sequence Similarities	Belongs to the LAMP family.
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

BACKGROUND

Introduction	Presents carbohydrate ligands to selectins. Also implicated in tumor cell metastasis.
Keywords	LAMP1; CD107 antigen like family member A; CD107 antigen-like family member A; CD107a; CD107a antigen; LAMP 1; LAMP-1; LAMPA; LGP120; IgpA; Lysosomal membrane glycoprotein 120KD; Lysosomal Associated Membrane Protein 1; Lysosome associated membrane glycop

GENE INFORMATION

Entrez Gene ID	3916
UniProt ID	P11279