



Human LLGL1 peptide (DAG-P2065)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Format	Lyophilised
Buffer	Lyophilized with 100% acetonitrile.
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles. Lyophilized with 100% acetonitrile.

BACKGROUND

Introduction	LLGL1 is a protein that is similar to a tumor suppressor in Drosophila. The protein is part of a cytoskeletal network and is associated with nonmuscle myosin II heavy chain and a kinase that specifically phosphorylates this protein at serine residues. The
Keywords	LLGL1 + LLGL2; DLG4; HGL; HUGL; HUGL 1; HUGL1; Human homolog to the D lgl gene protein; lethal giant larvae homolog 1; lethal giant larvae homolog 2; Lethal(2) giant larvae protein homolog 1; Lethal(2) giant larvae protein homolog 2; LGL2; LLGL

GENE INFORMATION

Entrez Gene ID	3996
UniProt ID	Q15334