



V5 tag peptide (DAG-P2075)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

BACKGROUND

Introduction	The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKIPNPLLGLDST), although it has also been used with a shorter 9 amin
Keywords	V5 tag; GKIPNPLLGLDST epitope tag; GKIPNPLLGLDST tag; Protein Rev; Regulator of expression of viral proteins; rev; V5 epitope tag