



## RFP peptide (DAG-P2077)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |   |
|--------------|---|
| Conjugate    | Unconjugated  |
| Format       | Liquid  |
| Preservative | None  |
| Storage      | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request. |

### BACKGROUND

|              |  |
|--------------|--|
| Introduction | Red fluorescent protein cloned from Discosoma coral (DsRed or drFP583) is an important tool in biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein (GFP; from Aequorea jellyfish). |
| Keywords     | RFP; DsRed; DsRed2; GFP like chromoprotein; red fluorescent protein; Red fluorescent protein drFP583; Red fluorescent protein eqFP611; Red fluorescent protein FP593; Red fluorescent protein zoanRFP  |