



ADV type 40 Hexon (DAG-P2192)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Adenovirus type 40 full length protein
Species	ADV
Purity	> 95% as determined by SDS-PAGE
Conjugate	Unconjugated
Format	Liquid
Size	100 µg
Buffer	10mM Bis Tris Propane pH7.2, 375 mM NaCl
Preservative	0.09% Sodium Azide
Storage	Store at -80°C

BACKGROUND

Introduction	Adenoviruses (members of the family Adenoviridae) are medium-sized (90–100 nm), nonenveloped (without an outer lipid bilayer) viruses with an icosahedral nucleocapsid containing a double stranded DNA genome. Their name derives from their initial isolation
Keywords	Astrovirus; Adv 40