



Toscana Virus Nucleocapsid (aa 1 - 180) (DAG-P2281)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Toscana Virus Nucleocapsid full length protein
Antigen Description	At least five proteins have been identified in Toscana virus-infected cells: nucleoprotein N, glycoproteins G1 and G2, a large protein (L) assumed to be a component of the polymerase, and two nonstructural proteins, NSm and NSs.
Species	Toscana Virus
Purity	>95% by SDS-PAGE .
Conjugate	Unconjugated
Applications	WB
Molecular Weight	34 kDa including tags
Format	Liquid
Buffer	pH: 8.00Constituents: 45% Phosphate Buffer, 48.5% Urea, 5.8% Sodium chloride, 0.1% C13E10
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 8.00Constituents: 45% Phosphate Buffer, 48.5% Urea, 5.8% Sodium chloride, 0.1% C13E10

BACKGROUND

Introduction	Toscana virus (TOSV) is an arbovirus (arthropod-borne virus) belonging to Bunyaviridae, a
---------------------	---

family of negative-stranded, enveloped RNA viruses. The virus can be transmitted to humans by the bite of an infected sandfly of the genus *Phlebotomus*. Toscana is no

Keywords

N; Nucleocapsid protein; Protein N; TOSV Nucleocapsid
