



HIV type 1 P311 (full length) (DAG-P2466)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	HIV-1 p31 full length protein
Species	HIV
Purity	>95% by SDS-PAGE .
Conjugate	Unconjugated
Applications	ELISA WB
Format	Liquid
Buffer	Information available upon request.
Preservative	None
Storage	Shipped at 4°C. Store at -20°C. Information available upon request.

BACKGROUND

Introduction	The human immunodeficiency virus (HIV) is a lentivirus (slowly replicating retrovirus) that causes the acquired immunodeficiency syndrome (AIDS), a condition in humans in which progressive failure of the immune system allows life-threatening opportunistic
Keywords	HIV1 p311