



Recombinant Dengue virus type 1 NS1 protein (DAG-P2497)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥95% as determined by SDS-PAGE
Conjugate	His
Applications	LFIA, CLIA
Format	Liquid
Size	200 µg, 1 mg
Buffer	PBS, 5% mannitol, 0.01% Tween 80, pH7.4
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Dengue virus; DENV; Nonstructural protein 1; NS1
----------	--