



Recombinant Dengue virus type 2 NS1 protein (DAG-P2500)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------|---|
| Purity | ≥95% as determined by SDS-PAGE |
| Conjugate | His |
| Applications | LFIA, CLIA |
| Format | Liquid |
| Size | 200 µg, 1 mg |
| Buffer | PBS, 5% mannitol, 0.01% Tween 80, pH7.4 |
| Storage | Store at -20°C to -80°C |

BACKGROUND

| | |
|----------|--|
| Keywords | Dengue virus; DENV; Nonstructural protein 1; NS1 |
|----------|--|