



HAV (aa 1392 - 1521) (DAG-P2692)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Hepatitis A Virus protein fragment
Species	HAV
Purity	> 90 % by SDS-PAGE.
Conjugate	Unconjugated
Applications	WB ELISA SDS-PAGE
Bio-activity	Reacts strongly with human Hepatitis A Virus positive serum.
Format	Liquid
Buffer	Preservative: None Constituents: 50% Glycerol, 0.1% SDS, 10mM Tris HCl, pH 9.6
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles. Preservative: None Constituents: 50% Glycerol, 0.1% SDS, 10mM Tris HCl, pH 9.6

BACKGROUND

Introduction	Hepatitis A (formerly known as infectious hepatitis) is an acute infectious disease of the liver caused by the hepatitis A virus (HAV). Many cases have little or no symptoms especially in the young. The time between infection and symptoms, in those who de
Keywords	HAV; Hepatitis A; Hepatitis A Virus