



## HAV (aa 1121 - 1234) (DAG-P2808)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Active Hepatitis A Virus protein fragment
<b>Species</b>	HAV
<b>Purity</b>	> 90 % by SDS-PAGE. This antigen was purified by a proprietary chromatographic technique.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB ELISA
<b>Bio-activity</b>	Active Hepatitis A Virus protein fragment is immunoreactive with sera from HAV-infected individuals.
<b>Format</b>	Liquid
<b>Buffer</b>	Preservative: None Constituents: 10mM CBB, 50% Glycerol, 0.1% SDS, pH 9.6
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze/thaw cycles.

### BACKGROUND

<b>Introduction</b>	Hepatitis A (formerly known as infectious hepatitis) is an acute infectious disease of the liver caused by the hepatitis A virus (HAV). Many cases have little or no symptoms especially in the young. The time between infection and symptoms, in those who de
<b>Keywords</b>	HAV; Hepatitis A; Hepatitis A Virus