



Recombinant Hantavirus Nucleocapsid (Hantavirus Dobrava, strain Slovenia) (DAG-P2859)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	55 kDa
Format	Liquid or Lyophilized powder
Size	100 µg, 1 mg
Buffer	PBS-based buffer
Storage	Store at 2-8°C

BACKGROUND

Keywords	Hantavirus; Nucleocapsid
-----------------	--------------------------