



Recombinant Hantavirus Nucleocapsid (Hantavirus Seoul, strain 80-39) (DAG-P2861)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity > 95% as determined by SDS-PAGE

Conjugate Unconjugated

Applications ELISA, WB

Molecular Weight 55 kDa

Format Liquid or Lyophilized powder

Size 100 µg, 1 mg

Buffer PBS-based buffer

Storage Store at 2-8°C

BACKGROUND

Keywords Hantavirus; Nucleocapsid