



Retroviruses Viral Protein U (DAG-P2931)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Viral Protein U protein fragment
Species	Retroviruses
Purity	> 95 % by SDS-PAGE.
Conjugate	Unconjugated
Applications	SDS-PAGE
Format	Liquid
Buffer	Preservative: None Constituents: HEPES, pH 7.2
Preservative	None
Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Preservative: None Constituents: HEPES, pH 7.2

BACKGROUND

Introduction	Retroviruses have several characteristic structural and catalytic proteins. One such auxiliary protein is a viral protein U (Vpu) which enhances virion release from human cells and is also involved in the degradation of CD4, the cellular surface receptor
Keywords	U ORF protein; Vpu