



Human CA 50 (DAG-T1007)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Colorectal Cancer, Gastrointestinal Malignancies
Species	Human
Purity	Control Grade
Conjugate	Unconjugated
Format	Liquid
Size	1 kU
Buffer	Liquid in 0.1M PBS buffer, pH 7.4
Preservative	0.05% Sodium Azide
Storage	-20°C

BACKGROUND

Introduction	CA50 is a carcinoma associated carbohydrate antigen.
Keywords	Antigen CA50; CA 50; CA 50 antigen; CA50 antigen; Carbohydrate antigen CA50; Tumour marker CA50; CA50