



Human Laminin (DAG-T1015)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	SDS PAGE: >90%
Conjugate	Unconjugated
Cas.No	114956-81-9
Molecular Weight	In electrophoresis a major band is observed at 60 kDa and a minor band at 170-190 kDa.
Format	Liquid
Concentration	0.5 mg/ml
Preservative	None
Storage	-20°C

BACKGROUND

Introduction	Laminins are major proteins in the basal lamina (formerly improperly called "basement membrane"), a protein network foundation for most cells and organs. The laminins are an important and biologically active part of the basal lamina, influencing cell differentiation, migration, adhesion as well as phenotype and survival.
Keywords	Laminin