



# Recombinant HBV Surface Antigen (subtype adw2, mutation G-145-R) (DAG-T1037)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	N/A
<b>Molecular Weight</b>	25 kDa (monomer)
<b>Format</b>	Liquid
<b>Size</b>	100 µg, 1 mg
<b>Buffer</b>	50 mM Tris-HCl, pH8.0, 200 mM NaCl
<b>Preservative</b>	None
<b>Storage</b>	Store at -80°C

## BACKGROUND

<b>Keywords</b>	HBsAg; HBV surface protein
-----------------	----------------------------