



Recombinant HBV Surface Antigen (subtype ayw1, mutation G-145-R) (DAG-T1049)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	N/A
Molecular Weight	25 kDa (monomer)
Format	Liquid
Size	100 µg, 1 mg
Buffer	50 mM Tris-HCl, pH8.0, 200 mM NaCl
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Dry ice

BACKGROUND

Keywords	HBsAg; HBV surface protein
----------	----------------------------