



CMV Pp38 (aa 117 - 373) (DAG-T1119)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	E.coli derived recombinant. The protein contains the CMV Pp38 (UL80a) immunodominant region (117-373aa).
Species	CMV
Conjugate	Unconjugated
Preservative	None
Storage	Store at 2–8 °C

BACKGROUND

Introduction	Cytomegalovirus (from the Greek cyto-, "cell", and -megalo-, "large") is a viral genus of the viral family known as Herpesviridae or herpesviruses. It is typically abbreviated as CMV.
Keywords	Cytomegalovirus Pp38; CMV Pp38; Cytomegalovirus; CMV