



Recombinant Cytomegalovirus Mosaic[His] (DAG-T1120)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|-------------------------|---|
| Product Overview | E.coli derived recombinant. The artificial mosaic protein contains the CMV gB immunodominant regions. |
| Species | CMV |
| Conjugate | Unconjugated |
| Preservative | None |
| Storage | Store at 2–8 °C |

BACKGROUND

| | |
|---------------------|---|
| Introduction | Cytomegalovirus (from the Greek cyto-, "cell", and -megalo-, "large") is a viral genus of the viral family known as Herpesviridae or herpesviruses. It is typically abbreviated as CMV. |
| Keywords | Cytomegalovirus; CMV; CytoMegalo mosaic; CMV mosaic |