



SARS Nucleocapsid (DAG-T1127)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	E.coli derived recombinant. The protein contains the nucleocapsid protein immunodominant regions.
Species	SARS
Conjugate	Unconjugated
Preservative	None
Storage	Store at 2–8 °C

BACKGROUND

Introduction	Severe acute respiratory syndrome (SARS) is a viral respiratory disease of zoonotic origin caused by the SARS coronavirus (SARS-CoV). Between November 2002 and July 2003, an outbreak of SARS in southern China caused an eventual 8,273 cases and 775 deaths reported in multiple countries with the majority of cases in Hong Kong (9.6% fatality rate) according to the World Health Organization (WHO). Within weeks, SARS spread from Hong Kong to infect individuals in 37 countries in early 2003.
Keywords	SARS nucleocapsid (N)