



La (SSB) Antigen (DAG-T1217)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Purified from bovine thymus. After coating onto ELISA plates the product will bind autoantibodies to La (SSB).
Purity	The La autoantigen (45-50 kDa) is more than 90% pure, as assessed by SDS polyacrylamide gel electrophoresis.
Concentration	0.1-1.0 mg/ml
Size	0.20 mg
Preservative	None
Storage	The product is stabilised with 20% glycerol and 0.1% Micr-O-protect TM. Store at -20°C or below (long term) or at +4°C (short term). Avoid repeated freezing and thawing. Mix thoroughly before use.

BACKGROUND

Introduction	Sjögren syndrome type B antigen (SS-B) also known as Lupus La protein is a protein that in humans is encoded by the SSB gene.
Keywords	SSB;Single-strand binding protein;La (SSB) Antigen;La (SSB);Antigen;La Antigen;SSB Antigen

GENE INFORMATION

Entrez Gene ID	6741
UniProt ID	P05455