



HIV Glycoprotein 41 [Gold] (DAG-T1424)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	rHIV gp41 antigen was conjugated to a 40 nm colloidal gold particle at 10 micrograms per mL
Species	HIV
Conjugate	Gold
Size	1 ml, 10 ml, 100 ml
Buffer	Bovine Serum Albumin/Phosphate, pH 7.4
Preservative	0.1% Sodium Azide
Storage	4°C

BACKGROUND

Introduction	HIV-1 is the most common and pathogenic strain of the virus. Scientists divide HIV-1 into a Group M (major) and two or more minor groups.
Keywords	HIV-1 gp41; HIV1 gp41 (HRP); Env; Env polyprotein; Envelope glycoprotein gp160; Envelope Protein gp120; Glycoprotein 41; GP120; gp41; Human Immunodeficiency Virus 1; Human Immunodeficiency Virus Type 1 gp41; TM; Transmembrane protein; Retroviridae; Lentiv