



Gamithromycin [KLH] (DAG-WT016K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Gamithromycin is a new type of semi-synthetic macrolide antibiotics for veterinary use. It is one of the representative drugs of the second generation of macrolide antibiotics. It inhibits bacterial RNA affects protein synthesis and exerts antibacterial and Sterilizing effect. Gamithromycin has high antibacterial activity against Mannella hemolyticus, Pasteurella multocida and Mycoplasma bovis, which cause bovine respiratory diseases.
--------------	--