



Acid Orange [BSA] (DAG-WT020)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Acid Orange conjugated with BSA
Specificity	This product is ELISA-paired with CABT-Z529M
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Acid orange is an organic compound and an azo dye. It is one of the first water soluble dyes to be commercialized and also be used as a pH indicator. It is produced by azo coupling of β -naphthol and diazonium derivative of sulfanilic acid.
--------------	--