



Penconazole [KLH] (DAG-WT024K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	1-[2-(2,4-dichlorophenyl)pentyl]1,2,4-triazole is a member of the class of triazoles that is 1,2,4-triazole substituted at position 1 by a 2-(2,4-dichlorophenyl)pentyl group. It is a dichlorobenzene and a member of triazoles.
--------------	---