



## Nitrofen [OVA] (DAG-WT047)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Product Overview	Nitrofen conjugated with OVA
Specificity	This product is ELISA-paired with CABT-Z543M
Species	N/A
Conjugate	OVA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Introduction	Nitrofen is an organic molecular entity. It has a role as an EC 1.3.3.4 (protoporphyrinogen oxidase) inhibitor and a herbicide.
--------------	---