



Eugenol [OVA] (DAG-WT057)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Eugenol conjugated with OVA
Specificity	This product is ELISA-paired with CABT-Z548M
Species	N/A
Conjugate	OVA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Eugenol is a phenylpropanoid formally derived from guaiacol with an allyl chain substituted para to the hydroxy group. It has a role as an allergen, a plant metabolite, a human blood serum metabolite and a sensitiser. It is a phenylpropanoid, a monomethoxybenzene and a member of phenols. It derives from guaiacol.
---------------------	--