



Ranitidine [BSA] (DAG-WT060)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Ranitidine conjugated with BSA
Specificity	This product is ELISA-paired with CABT-Z549M
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Ranitidine is a member of the class of furans used to treat peptic ulcer disease (PUD) and gastroesophageal reflux disease. It has a role as an anti-ulcer drug, a H2-receptor antagonist, an environmental contaminant, a xenobiotic and a drug allergen. It is a member of furans, a tertiary amino compound, a C-nitro compound and an organic sulfide.
---------------------	--