



2,4-Toluene diisocyanate [BSA] (DAG-WT064)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	2,4-Toluene diisocyanate conjugated with BSA
Specificity	This product is ELISA-paired with CABT-Z551M
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Toluene diisocyanate exists in two isomeric forms (2,4-toluene diisocyanate and 2,6-toluene diisocyanate) which have similar properties and effects. Toluene diisocyanate is produced commercially as an 80:20 (2,4-toluene diisocyanate:2,6-toluene diisocyanate) mixture of the two isomers.
---------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------