



Hyoscyamine [BSA] (DAG-WT068)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Hyoscyamine conjugated with BSA
Specificity	This product is ELISA-paired with CABT-Z553M
Species	N/A
Conjugate	BSA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Hyoscyamine is a chemical compound, a tropane alkaloid it is the levo-isomer to atropine. It is a secondary metabolite of some plants, particularly henbane (<i>Hyoscamus niger</i> .) Hyoscyamine is used to provide symptomatic relief to various gastrointestinal disorders including spasms, peptic ulcers, irritable bowel syndrome, pancreatitis, colic and cystitis.
--------------	--