



Triclosan [OVA] (DAG-WT073)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Triclosan conjugated with OVA
Specificity	This product is ELISA-paired with CABT-Z556M
Species	N/A
Conjugate	OVA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Triclosan is an aromatic ether that is phenol which is substituted at C-5 by a chloro group and at C-2 by a 2,4-dichlorophenoxy group. It has a role as an antibacterial agent, a fungicide, an EC 1.5.1.3 (dihydrofolate reductase) inhibitor, an antimalarial, an EC 1.3.1.9 [enoyl-[acyl-carrier-protein] reductase (NADH)] inhibitor, a xenobiotic, a persistent organic pollutant and a drug allergen.
---------------------	---