



DNP [OVA] (DAG-WT089)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	DNP conjugated with OVA
Species	N/A
Conjugate	OVA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	5 mg, 10 mg, 100 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Store at <4°C. Protect from light

BACKGROUND

Introduction	A toxic dye, chemically related to trinitrophenol (picric acid), used in biochemical studies of oxidative processes where it uncouples oxidative phosphorylation. It is also used as a metabolic stimulant.
Keywords	2,4-Dinitrophenol; DNP; Dinitrophenols