



Recombinant Human MMP-3 [His] (DAG-WT101)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	53.3 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM Tris-HCl, 150mM NaCl, 0.05% Brij35, 10% Glycerol, pH 7.5
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Stromelysin-1 also known as matrix metalloproteinase-3 is an enzyme that in humans is encoded by the MMP3 gene. The MMP3 gene is part of a cluster of MMP genes which localize to chromosome 11q22.3. MMP-3 has an estimated molecular weight of 54 kDa.
Keywords	Matrix metalloproteinase-3; MMP

GENE INFORMATION

