



Gastrin-17 Antigen (DAG-WT102)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	>95% by SDS-PAGE
Conjugate	Unconjugated
Applications	Control
Reconstitution	Reconstitute with 10 mM PBS, pH 7.4
Format	Lyophilized
Size	1 mg
Preservative	None
Storage	Store at $\leq -20^{\circ}\text{C}$. Avoid repeated freezing and thawing.
Ship	Cold Packs

BACKGROUND

Introduction	Little gastrin I is a form of gastrin commonly called as gastrin-17. This is a protein hormone, secreted by the intestine. Gastrin II has identical amino acid composition to Gastrin I, the only difference is that the single tyrosine residue is sulfated in Gastrin II.
Keywords	Little gastrin I; Gastrin-17; G-17