



alpha-1-Microglobulin Antigen (alpha-1-MG Ag) (DAG-WT104)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Control/Calibrators
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	20mM PBS, 350mM NaCl, 500mM Imidazole, 10% Glycerol, pH7.2
Preservative	0.05% Proclin300
Storage	2-8 °C (up to 1 Year). Avoid freezing

BACKGROUND

Introduction	Alpha-1-microglobulin is a microglobulin, a small globular protein. It is found in all vertebrates, including humans, and is distributed in blood plasma and extravascular tissues of all organs. It is synthesized in most cells of the body, but mainly in the liver from a gene that codes for the alpha-1-microglobulin/bikunin precursor.
--------------	--