



Phenolphthalein [KLH] (DAG-WT1068K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Phenolphthalein is an organic compound which is used medicinally as a cathartic. It has a chemical formula $C_{20}H_{14}O_4$. It is written as "phph" or HIn". It is used in acid-base titration as an indicator. As an indicator it turns pink to red in alkaline and is colourless in acid solutions.
--------------	--