



Cholyglycine-BSA Antigen (CG-BSA Ag) (DAG-WT107)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	Control/Calibrators
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	20mM PBS
Preservative	0.05% Proclin300
Storage	2-8 °C (up to 1 Month), below -20 °C (up to 2 Years). Avoid repeat freeze thaw cycles.

BACKGROUND

Introduction	Glycocholic acid is a bile acid glycine conjugate having cholic acid as the bile acid component. It has a role as a human metabolite. It derives from a cholic acid and a glycochenodeoxycholic acid. It is a conjugate acid of a glycocholate.
--------------	---