



## Dimethyl p-nitrophenylphosphate [BSA] (DAG-WT1076B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                                                                       |
|---------------|---------------------------------------------------------------------------------------|
| Conjugate     | BSA                                                                                   |
| Applications  | ELISA, LFIA                                                                           |
| Format        | Liquid                                                                                |
| Concentration | Batch dependent - please inquire should you have specific requirements                |
| Size          | 1 mg                                                                                  |
| Buffer        | 10mM PBS, pH7.4                                                                       |
| Preservative  | None                                                                                  |
| Storage       | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| Ship          | Wet ice                                                                               |