



Etoxazole [BSA] (DAG-WT1077B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Ettoxazole is a narrow spectrum systemic acaricide used to combat spider mites. It targets a variety of mites in the egg, larvae and nymph stages however not the adult stage. It also exhibits insecticidal activity towards aphids, the green rice leafhopper and diamondback moth.
--------------	---