



Recombinant Infectious Bronchitis Virus N Protein (DAG-WT1121)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90 % as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	50mM Tris, pH 8.0, 300mM NaCl, 2mM DTT, 2mM EDTA, 0.1% Triton X-100
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Infectious bronchitis virus (IBV) is a highly infectious avian pathogen which affects the respiratory tract, gut, kidney and reproductive systems. It belongs to the Coronavirus family, genus Gammacoronavirus which has a single stranded RNA genome, surrounded by a nucleocapsid and envelope.
Keywords	Infectious bronchitis virus; IBV; Nucleocapsid