



Recombinant PRRSV GP5 Protein [His] (DAG-WT1143)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, WB
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS buffer
Preservative	None
Storage	Store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction Porcine reproductive and respiratory syndrome virus (PRRSV) is one of the most significant etiological agents in the swine industry worldwide. It has been reported that PRRSV infection can modulate host immune responses, and innate immune evasion is thought to play a vital role in PRRSV pathogenesis.

Keywords Porcine reproductive and respiratory syndrome virus; PRRSV; GP5

GENE INFORMATION

