



Recombinant HIV Reverse Transcriptase (DAG-WT1145)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	N/A
Bio-activity	≥10,000 units/mL ≥5,000 units/mg protein
Format	Liquid
Size	5000 U
Buffer	10 mM potassium phosphate, 1 mM DTT, 20% glycerol, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	HIV-1 reverse transcriptase is an RNA-dependent DNA polymerase of HIV-1(AIDS virus), subtype B origin. It also has RNaseH activity and is an enzyme indispensable for reproduction of AIDS virus.c
Keywords	Human immunodeficiency virus; HIV; Reverse transcriptase