



Apolipoprotein B Antigen (ApoB Ag) (DAG-WT115)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunogen/Calibrators/Control
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	20mM PBS, pH7.2
Preservative	0.05% Proclin300
Storage	2-8 °C (up to 1 Year). Avoid freezing

BACKGROUND

Introduction	Apolipoprotein B is the primary apolipoprotein of chylomicrons, VLDL, Lp(a), IDL, and LDL particles (LDL - known commonly by the misnomer "bad cholesterol" when in reference to both heart disease and vascular disease in general), which is responsible for carrying fat molecules (lipids), including cholesterol, around the body to all cells within all tissues.
--------------	---