



# Pneumonia Virus of Mice (PVM) Positive Control (DAG-WT1153)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Mouse
Conjugate	N/A
Applications	Control
Format	Liquid
Size	1 mL
Storage	Store at 4°C
Ship	Wet ice

## BACKGROUND

Introduction	Pneumonia Virus of Mice (PVM) is related to the human and bovine respiratory syncytial virus (RSV) pathogens, and has been used to study respiratory virus replication and the ensuing inflammatory response as a component of a natural host—pathogen relationship.
Keywords	Pneumonia Virus of Mice; PVM